


- **Still using ISPyB GUI in production (Version 5.13.1)**
- **Discussions that were ongoing internally at SOLEIL ended up by the decision to switch directly to EXI2**
- **Discussions with SOLEIL direction allowed us to substract the work for EXI2 installation in production : to be done in July**
- **We are going to set up an open source Single Sign-On – SSO (KeyCloak):**
  - **EXI / ISPyB**
  - **The User Office Tools (SUN set)**
  - **...**
- **Used for MX and SAXS (a priori no needs for other techniques)**

## Issue #52 : Add information about HPLC column in the column "Run" of the HPLC results panel

Add three new parameters, that could fit in the column Run of both buffer and sample lines, just below the "#": "Column name" (ex: Biosec300) "Column Volume" (ex : 5 ml) "Elution Flux" (ex: 0.3 ml/min)



## Issue #53 : Adding the UV absorption profile in addition to the already present Rg/I0/Mass profiles of an HPLC run

